

Steps for conducting intestinal ultrasound

Federica Furfaro

Assistant Professor (Tenure Track) Vita-Salute University San Raffaele Gastroenterology and Endoscopy Milan, Italy



Kuwait City, Kuwait, November 21st – 22nd,2025





Disclosure

I have received:

- •Speaker fees from Janssen, AbbVie, Pfizer, Galapagos, Tillotts, Alfa Sigma
- •Advisory board/consultancy fees from Janssen, AbbVie, Pfizer, Galapagos





Prepare the machine and yourself



- Know the indication, clinical setting and history of the patient
- Clean hands/probe (NO alcoholic solutions)
- On the right side of the patient
- Cord does not belong on the floor or around your neck
- Probe in right hand
- Left hand on buttons



Menti



Join at menti.com | use code 7589 3778


Is fasting a preparation requirement for the patient?




Mentimeter


Mod 1 kuwait_Steps for ...

Choose a slide to present



Select a different Menti





Why are landmarks important?



Always required

Not required

Preferrred for motility assessment

→ Start quiz





Preparation (patient)



Fasting


- not necessary
- ->4 h -> reduce amount of food and air in the small bowel
- ->6 h → recommended before assessing splanchnic blood flow and GI motility





Prepare the machine – Insert patient's details

- Check patient ID / details
- Save examination when you finish







Ultrasound machine



Monitor

Touchscreen

Probes

Control panel

Gel








Ultrasound probes





Frequency





Properties of the different probes

Low frequency: 1-6 MHz,

1-6 MHz, overview, more depth, less resolution


curved array for a larger field, dependent pathology



- For a complete examination of the bowel both a low and high resolution probe are needed, LoE 5, GoR C, Strong consensus 13/13
- A probe with a frequency above 5 MHz should be used when measuring wall thickness, LoE 4, GoR B, Strong consensus 13/13

Nylund K et al. EFSUMB Recommendations and Guidelines for Gastrointestinal Ultrasound, Ultraschall Med 2017

High frequency: > 5MHz (5-10 MHz), only surface accessible (8-10 cm depth



penetration), high resolution, detailed examination (discrimination of stratification of the bowel wall)







Ultrasound technique










Prepare the patient












- Expose abdominal area adequately
- Water-based ultrasound gel
- Hold the probe correctly








Ultrasound technique






In a segment-wise manner, **all portions** of the bowel are scanned in long and short axes, and images are simultaneously captured to ensure **a continuous assessment**






Scan planes



Imaging




Intestinal ultrasound







Scan planes

Longitudinal





Cross-section









Probe orientation



caudal








Menti


Mod1kuwait_Steps for ... [♂






Choose a slide to present








Why are landmarks important?







Turning the probe
























Keep your region of interest in the center











Keep your region of interest in the center











Keep your region of interest in the center











Examination technique

- Probe movement
 - → set and tilt, not chaotic movement
- Graded compression method
 - → visualization of deep areas using HF probe











What to look for

- Recognize the different parts of the bowel and surrounding structures
 - → Landmarks!!


Normal vs. abnormal bowel wall





Localization of the key anatomy (landmarks)

The most important landmarks: iliac vessels and the psoas muscle












Anatomical landmarks in the mid upper andomen - Finding the Transverse







Waiting for players...















Normal bowel wall

5 distinct sonographic layers










Normal bowel wall in ultrasonography



Some practical points:

- Interface echo of serosa mixed in with interface from peritoneum
- Posterior bowel wall often not possible to see due to air in the lumen → Measurements in anterior wall!
- Measure from mucosal hypoechoic layer to muscularis hypoechoic layer
- Avoid measuring at points of haustrations and mucosal folds (overestimation!)
- Keep probe perpendicular to GI wall





Bowel wall thickness (BWT)



^{*}Measure at right angles to bowel wall layers.




Parameters




- ✓ BWT
- ✓ Stratification
- √ Vascolarization
- ✓ Margins
- ✓ Austra
- ✓ Motility



- ✓ Mesentery
- ✓ Lymph nodes
- ✓ Free fluid












Examination technique – Systematic


RECOMMENDATIONS:

 The scanning of the intestines must involve a systematic approach, LoE 5, GoR C. Strong consensus 12/12 Rectum → Sigmoid → Decending → Transverse → Ascending → Ceacum


- → Terminal Ileum
- → Small bowel

Colon





Rectum







Sigmoid Colon









Sigmoid Colon












Descending colon









Transverse colon










Ascending colon









Terminal ileum








Examination technique Systematic approach

Small bowel

"mowing the lawn" technique











Differentiation between small and large bowel on IUS









Take home messages

- ✓ Preparation: none
- ✓ Probe 5-8 MHz
- ✓ orientation

- ✓ know what to look for
- ✓ anatomical landmarks
- ✓ ultrasound scan plans
- ✓ systematic approach



